

TSV™ BellaTek® Encode® Impression System



TSV BellaTek Encode Healing Abutments

Surgical Placement

The **TSV BellaTek Encode Healing Abutment** is a two-piece healing abutment designed to facilitate gingival tissue healing and consists of an abutment and retaining screw that are assembled together.

TSV BellaTek Encode Healing Abutments have the Encode Coding scheme on the occlusal surface and a hex connection at the base of the healing abutment that engages the hex with the implant for orientation and anti-rotation.

Note: Please ensure that the instructions for use are made available to the restorative clinician and use the patient chart labels located in the package of the BellaTek Encode Healing Abutment.



For instructions for use: ifu.biomet3i.com



Following implant placement, select a TSV BellaTek Encode Healing Abutment with the appropriate restorative platform diameter, emergence profile and collar height.



Use wide TSV BellaTek Encode Healing Abutment(s) to ensure the final abutment will seat easily without excessive blanching of the soft tissue. Bone profiling may be required to fully seat the healing abutment, particularly in cases with sub-crestal implant placement.



When possible use tall TSV BellaTek Encode Healing Abutments. The height of the abutment collar, not including the domed occlusal portion, should extend 1-2 mm above the soft tissue (1 mm minimum) on all sides to facilitate proper traditional or digital impression and subsequent scanning.




Pick up the TSV BellaTek Encode Healing abutment(s) using a hexed GemLock® driver (HXLGR1.25 or HXLGR1.25) and place the healing abutment on the implant and hand-tighten.

TSV™ BellaTek® Encode® Impression System



TSV BellaTek Encode Healing Abutments



Following placement, confirm that the TSV BellaTek Encode Healing Abutment is fully seated on the implant with a radiograph.



Once fully seated, torque the healing abutment screw to 20 Ncm using a restorative torque wrench (TWR) and torque wrench hex driver (TW1.25 or TW1.25L). Where applicable, suture the soft tissue around the TSV BellaTek Encode Healing Abutment.

Ordering Information



Collar Height	3.5 mm (D)		4.5 mm (D)			5.7 mm (D)
	3.8 mm (EP)	5.0 mm (EP)	5.0 mm (EP)	5.6 mm (EP)	6.0 mm (EP)	6.8 mm (EP)
3 mm	TEHA3383	TEHA3503	TEHA4503	TEHA4563	TEHA4603	TEHA5683
5 mm	TEHA3385	TEHA3505	TEHA4505	TEHA4565	TEHA4605	TEHA5685
7 mm	TEHA3387	--	--	TEHA4567	--	--

D: Restorative Platform Diameter • EP: Emergence Profile

Contact us at 1-800-342-5454 or visit zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZBINST0025 REV B 01/20 ©2020 Zimmer Biomet. All rights reserved.

