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EN ISO 13485:2016 + A11:2021 – Medical Devices – Quality management systems – Requirement for regulatory purposes

EN ISO 10993-1:2020 – Biological evaluation of medical Devices – Part 1: Evaluation and testing within a risk management process

EN ISO 14971:2019 + A11:2021 – Medical Devices – Application of risk management to medical devices

EN ISO 11137-1:2015 + A2:2019 – Sterilization of health care products – Radiation – Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices

EN ISO 11137-2:2015+A1:2023 – Sterilization of health care products – Radiation – Part 2: Establishing the sterilization dose

EN 556-1:2001 – Sterilization of medical devices – Requirements for medical devices to be designated as “STERILE” – Part 1: Requirements for terminally sterilization medical devices

EN ISO 20417:2021 – Medical devices - Information to be supplied by the manufacturer

EN ISO 11607-1:2020 + A11:2022 – Packaging for Terminally Sterilized Medical Devices – Part 1: Requirements for Materials, Sterile Barrier Systems and Packaging Systems

EN ISO 11607-2:2020 + A1:2023 – Packaging for Terminally Sterilized Medical Devices - Part 2: Validation Requirements For Forming, Sealing and Assembly Processes (ISO 11607-2:2006)

EN ISO 15223-1:2021+A1:2025 – Medical Devices – Symbols to be Used with Information to be Supplied by the Manufacturer - General Requirements

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10. Literatura

Gallucci GO, Hamilton A, Zhou W, Buser D, Chen S. Implant placement and loading protocols in partially edentulous patients: A systematic review. Clin Oral Implants Res. 2018; 29:106-34.

Ključna klinička ispitivanja
Abdelaal M, Fayyad A, Sheta N, AbdelNabi N, ELFar MM. Oral Health-Related Quality of Life in Single Implant Mandibular Overdenture Retained by CM LOC versus Ball Attachment: A Randomized Controlled Trial. Open Access Maced J Med Sci. 2019;7(21):3642-3646.

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Abdelaal, N.A., 2017. Randomized study on effect of conventional versus Schlosser's fabrication techniques for complete over dentures supported by mandibular implants. International Dental & Medical Journal of Advanced Research, 2017;3(1):1-5.
Agnini A, Salama MA, Agnini AM, Salama H, Stappert CF, Romeo D. RevitaliZe Patient Solutions: preliminary results from a single cohort prospective study using Screw-Vent TSVT implants. Dentalxp.com. Accessed December 26, 2022.
Akin R, Chapple AG. Clinical Advantages of Immediate Posterior Implants with Custom Healing Abutments: Up to 8-Year Follow-Up of 115 Cases. J Oral Maxillofac Surg. 2022;80(12):1952-1965.
Al-Magaleh WR, Swelem AA, Radi IAW. The effect of 2 versus 4 implants on implant stability in mandibular overdentures: A randomized controlled trial. J Prosthet Dent. 2017;118(6):725-731.
Arghami A, Simmons D, St Germain J, Maney P. Immediate and early loading of hydrothermally treated, hydroxyapatite-coated dental implants: a 7-year prospective randomized clinical study. Int J Implant Dent. 2021;7(1):21. Published 2021 Mar 10.
Bianconi S, Bozzoli P, Del Fabbro M. Treatment of Postextraction Sites With Allograft-Stabilized Dental Implants: A Clinical Case Series. Implant dentistry.2017;26(1):37-45.
Bianconi S, Wang HL, Testori T, Fontanella F, Del Fabbro M. Bone modifications around porous trabecular implants inserted with or without primary stability 2 months after tooth extraction: A 3-year controlled trial. Int J Oral Implantol (Berl). 2020;13(3):241-252.
Brauner E, Di Carlo S, Ciolfi A, et al. Use of Porous Implants for the Prosthetic Rehabilitation of Fibula Free Flap Reconstructed Patients. J Craniofac Surg. 2019;30(4):1163-1169.
Brauner E, Guarino G, Jamshir S, et al. Evaluation of Highly Porous Dental Implants in Postablative Oral and Maxillofacial Cancer Patients: A Prospective Pilot Clinical Case Series Report. Implant dentistry. 2015;24(5):631-637.
Brauner E, Jamshir S, Di Carlo S, Pagnoni M, Guarino G, Pompa G. Immediate implant loading: A comparison of Trabecular Metal and Tapered Screw-Vent dental implants. OHDM. 2015;14(2):1-6.
Cannizzaro G, Felice P, Minciarelli AF, Leone M, Viola P, Esposito M. Early implant loading in the atrophic posterior maxilla: 1-stage lateral versus crestal sinus lift and 8 mm hydroxyapatite-coated implants. A 5-year randomised controlled trial. Eur J Oral Implantol. 2013;6(1):13-25.
De Francesco M, Gobbato EA, Noce D, Cavallari F, Fioretti A. Clinical and radiographic evaluation of single tantalum dental implants: a prospective pilot clinical study. Oral Implantol (Rome). 2017;9(Suppl 1/2016 to N 4/2016):38-44.
Dimaira M. Immediate Placement of Trabecular Implants in Sites of Failed Implants. The International journal of oral & maxillofacial implants. 2019;34(5):e77-e83.
Dursun E, Tulunoglu I, Ozbek SM, et al. The influence of platform switching on clinical, laboratory, and image-based measures: a prospective clinical study. Clin Implant Dent Relat Res. 2014;16(6):936-946.

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Edelmann AR, Patel D, Allen RK, Gibson CJ, Best AM, Bencharit S. Retrospective analysis of porous tantalum trabecular metal-enhanced titanium dental implants. J Prosthet Dent. 2019;121(3):404-410.
El Chaar E, Castano A. A Retrospective Survival Study of Trabecular Tantalum Implants Immediately Placed in Posterior Extraction Sockets Using a Flapless Technique. The Journal of oral implantology. 2017;43(2):114-124.
El-Chaar ES. Immediate placement and provisionalization of implant-supported, single-tooth restorations: a retrospective study. Int J Periodontics Restorative Dent. 2011;31(4):409-419.
Farahat A, Nawar N. Comparative study of bone height changes around immediately loaded porous tantalum parallel sided Trabecular & Screw shaped implants retaining mandibular implant-overdentures using CBCT (RCT). Egypt Dent J. 2019;65(4):3571-3580.
Harel N, Piek D, Livne S, Palti A, Ormianer Z. A 10-year retrospective clinical evaluation of immediately loaded tapered maxillary implants. Int J Prosthodont. 2013;26(3):244-249.
Jang HW, Kang JK, Lee K, Lee YS, Park PK. A retrospective study on related factors affecting the survival rate of dental implants. J Adv Prosthodont. 2011;3(4):204-215.
Jung JH, Kim SY, Yi YJ, Lee BK, Kim YK. Hydroxyapatite-coated implant: Clinical prognosis assessment via a retrospective follow-up study for the average of 3 years. J Adv Prosthodont. 2018;10(2):85-92.
Keith JD Jr, Petrungaro P, Leonetti JA, et al. Clinical and histologic evaluation of a mineralized block allograft: results from the developmental period (2001-2004). Int J Periodontics Restorative Dent. 2006;26(4):321-327.
Khayat PG, Arnal HM, Tourbah BI, Sennerby L. Clinical outcome of dental implants placed with high insertion torques (up to 176 Ncm). Clin Implant Dent Relat Res. 2013;15(2):227-233.
Khayat PG, Milliez SN. Prospective clinical evaluation of 835 multithreaded tapered screw-vent implants: results after two years of functional loading. J Oral Implantol. 2007;33(4):225-231.
Kim JM, Sohn DS, Bae MS, Moon JW, Lee JH, Park IS. Flapless transcrestal sinus augmentation using hydrodynamic piezoelectric internal sinus elevation with autologous concentrated growth factors alone. Implant Dent. 2014;23(2):168-174.
Kim YK, Ahn KJ, Yun PY, et al. Effect of loading time on marginal bone loss around hydroxyapatite-coated implants. J Korean Assoc Oral Maxillofac Surg. 2013;39(4):161-167.
Koutouzis T, Huwais S, Hasan F, Trahan W, Waldrop T, Neiva R. Alveolar Ridge Expansion by Osseodensification-Mediated Plastic Deformation and Compaction Autografting: A Multicenter Retrospective Study. Implant Dent. 2019;28(4):349-355.
Lang LA, Edgin WA, Garcia LT, et al. Comparison of implant and provisional placement protocols in sinus-augmented bone: a preliminary report. Int J Oral Maxillofac Implants. 2015;30(3):648-656.
Lang LA, Turkyilmaz I, Edgin WA, Verrett R, Garcia LT. Immediate restoration of single tapered implants with nonoccluding provisional crowns: a 5-year clinical prospective study. Clin Implant Dent Relat Res. 2014;16(2):248-258.

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Lee CY, Rohrer MD, Prasad HS. Immediate loading of the grafted maxillary sinus using platelet rich plasma and autogenous bone: a preliminary study with histologic and histomorphometric analysis. <i>Implant Dent.</i> 2008;17(1):59-73.
Malmstrom H, Xiao J, Romanos G, Ren YF. Two-Year Success Rate of Implant-Retained Mandibular Overdentures by Novice General Dentistry Residents. <i>J Oral Implantol.</i> 2015;41(3):268-275.
Minichetti JC, D'Amore JC, Hong AY. Three-year analysis of Tapered Screw-Vent implants placed into extraction sockets grafted with mineralized bone allograft. <i>J Oral Implantol.</i> 2005;31(6):283-293.
Mostafa TM, El-Sheikh MM, Abd El-Fattah F. Implant-connected versus tooth-connected implant-supported partial dentures: 2-year clinical and radiographic comparative evaluation. <i>Int J Periodontics Restorative Dent.</i> 2015;35(3):335-343.
Nader N, Aboulhosn M, Berberi A, Manal C, Younes R. Marginal Bone Remodeling around healing Abutment vs Final Abutment Placement at Second Stage Implant Surgery: A 12-month Randomized Clinical Trial. <i>J Contemp Dent Pract.</i> 2016;17(1):7-15.
Nagi SE, Khan FR, Ali R. A 6-year Evaluation of 223 Tapered Dental Implants and associated prosthesis in 92 patients at a University Hospital. <i>J Pak Med Assoc.</i> 2016;66(Suppl 3)(10):S33-S35.
Ormianer Z, Palti A. Long-term clinical evaluation of tapered multi-threaded implants: results and influences of potential risk factors. <i>J Oral Implantol.</i> 2006;32(6):300-307.
Ormianer Z, Piek D, Livne S, et al. Retrospective clinical evaluation of tapered implants: 10-year follow-up of delayed and immediate placement of maxillary implants. <i>Implant Dent.</i> 2012;21(4):350-356.
Ormianer Z, Schirolli G. Maxillary single-tooth replacement utilizing a novel ceramic restorative system: results to 30 months. <i>J Oral Implantol.</i> 2006;32(4):190-199.
Peron C, Romanos G. Immediate Placement and Occlusal Loading of Single-Tooth Restorations on Partially Threaded, Titanium-Tantalum Combined Dental Implants: 1-Year Results. <i>The International journal of periodontics & restorative dentistry.</i> 2016;36(3):393-399.
Perret F, Romano F, Ferrarotti F, Aimetti M. Occlusive Titanium Barrier for Immediate Bone Augmentation of Severely Resorbed Alveolar Sockets with Secondary Soft Tissue Healing: A 2-Year Case Series. <i>Int J Periodontics Restorative Dent.</i> 2019;39(1):97-105.
Rahal RW, Shokry MM, Aboelsaad NS. Efficiency of Ridge Expansion Using Screw-Type Expanders with Simultaneous Trabecular Implant Placement in Narrow Anterior Maxilla (A Clinical and Radiographic Study) <i>Egyptian Dental Journal.</i> 2018;64(113):122.
Ranaan J, Bassir SH, Andrada L, et al. Clinical efficacy of the graft free slit-window sinus floor elevation procedure: A 2-year randomized controlled clinical trial. <i>Clin Oral Implants Res.</i> 2018;29(11):1107-1119.
Randall EF, Abou-Arraj RV, Geurs N, Griffin R, Reddy M, Geisinger M. The Effect of Dental Implant Collar Design on Crestal Bone Loss at 1 Year After Implant Placement. <i>Int J Periodontics Restorative Dent.</i> 2019;39(2):165-173.

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Schlee M, Pradies G, Mehmke WU, et al. Prospective, Multicenter Evaluation of Trabecular Metal-Enhanced Titanium Dental Implants Placed in Routine Dental Practices: 1-Year Interim Report From the Development Period (2010 to 2011). Clinical implant dentistry and related research. 2015;17(6):1141-1153.

Siddiqui AA, O'Neal R, Nummikoski P, et al. Immediate loading of single-tooth restorations: one-year prospective results. J Oral Implantol. 2008;34(4):208-218.

Simmons DE, Palaiologou A, Teitelbaum AG, Billiot S, Popat LJ, Maney P. Immediate and Early Loading of Hydrothermally Treated, Hydroxyapatite-Coated Dental Implants: 2-Year Results from a Prospective Clinical Study. Journal of Oral Implantology (2016) 42 (1): 17–25.

van der Schoor P, Schlee M, Wen H-B. Prospective Pilot Study of Immediately Provisionalized Restorations of Trabecular Metal-Enhanced Titanium Dental Implants: A 5-Year Follow-Up Report. Applied Sciences. 2022; 12(3):942.

Wen SC, Fu JH, Wang HL. Effect of Deproteinized Bovine Bone Mineral at Implant Dehiscence Defects Grafted by the Sandwich Bone Augmentation Technique. Int J Periodontics Restorative Dent. 2018;38(1):79-85.

Yoncheva E, Filtchev D. Comparison of the osseointegration of implants in smokers and non-smokers. Scr Sci Med Dent. 2020;6(1):13-17.

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ZimVie US Corp LLC,