

## Mobi-C® Cervical Disc Compared to M6-C™ Artificial Cervical Disc




# Technology Attributes

Table 1. Design & Biomechanics of Mobi-C and M6-C



**Mobi-C Cervical Disc**



**M6-C Artificial Cervical Disc**

Motion Technology: Mobile Bearing Core	Motion Technology: Deformable Core
<p><b>Technology Attributes</b></p> <ul style="list-style-type: none"> <li>• Mobile core adapts to the patient's natural movement through a 3-piece design intended to allow independent and coupled motion</li> <li>• Domed surface designed to articulate angularly with the superior endplate</li> <li>• Flat bottom designed to translate up to 1 mm and rotate on the inferior endplate</li> </ul>	<p><b>Technology Attributes</b></p> <ul style="list-style-type: none"> <li>• Incorporates an artificial nucleus (polycarbonate urethane) and a woven fiber annulus (polyethylene) designed to provide the same motion as a natural disc</li> </ul>
<p><b>Fixation Technology</b></p> <ul style="list-style-type: none"> <li>• Bone sparing lateral teeth designed to maintain integrity of endplates</li> <li>• Plasma sprayed titanium and hydroxyapatite coating and an anatomical dome on the superior endplate</li> </ul>	<p><b>Fixation Technology</b></p> <ul style="list-style-type: none"> <li>• 3 keels (2 mm high) on each endplate</li> <li>• Titanium plasma spray coating on the endplates</li> </ul>
<p><b>Range of Motion</b></p> <ul style="list-style-type: none"> <li>• Flexion/Extension: up to <math>\pm 10^\circ</math></li> <li>• Lateral bending: up to <math>\pm 10^\circ</math></li> <li>• Core rotation: up to <math>\pm 8^\circ</math></li> <li>• Core translation: <math>\pm 1</math> mm/axis</li> </ul>	<p><b>Range of Motion</b></p> <p>6 degrees of mobility</p> <ul style="list-style-type: none"> <li>• Flexion/Extension</li> <li>• Lateral bending</li> <li>• Axial rotation</li> <li>• Translation (A/P, lateral, vertical)</li> </ul>

## Clinical Analysis

Both Mobi-C and M6-C have independent studies with 2-year clinical data showing non-inferiority to ACDF at 1-level,<sup>1,2</sup> though there were notable differences between the studies in regards to patient inclusion criteria (Table 2) and reported outcomes (Tables 3 & 4).

Table 2. Notable Differences in the Inclusion/Exclusion Criteria of the 1-Level Mobi-C IDE Study and the 1-Level M6-C SSED Study<sup>1,2</sup>

	Mobi-C	M6-C
<b>Age</b>	18-75 years	18-75 years
<b>Disc Height</b>	$\geq 3$ mm	$\geq 4$ mm
<b>Severity of Facet Joint Arthritis</b>	Excluded radiographically confirmed severe facet joint disease or degeneration	Excluded symptomatic facet arthrosis
<b>DEXA Bone Mineral Density Score Cut-Off for Inclusion</b>	T-score $\geq -1.5$	T-score $> -2.5$
<b>Smoking Status</b>	Excluded $> 1$ pack/day	Did not exclude based on smoking status
<b>Diabetic Status</b>	Did not exclude diabetics	Excluded insulin-dependent diabetics
<b>VAS Inclusion Required</b>	Neck and/or arm pain VAS $\geq 30$ on 0-100 mm scale	Neck or arm pain VAS $\geq 4$ on 0-10 cm scale
<b>Prior Cervical Spine Surgery</b>	No prior spine surgery at the operative level, and no prior cervical fusion procedure at any level	No previous anterior cervical spine surgery, and no previous posterior cervical spine surgery that destabilizes the cervical spine
<b>NSAID Use</b>	Requested that patients refrain from prophylactic use of NSAIDs until 3 months postop	Allowed use of NSAIDs for 6 weeks postoperatively to improve HO

In the FDA IDE study of Mobi-C, 2-year follow-up showed Mobi-C to be non-inferior to ACDF at 1-level (Table 3).<sup>1</sup> In the 1-level arm of this prospective, randomized, controlled trial, Mobi-C had a 2-year success rate of 73.7%, compared to 65.3% for ACDF. Compared to ACDF, Mobi-C supported greater improvements in NDI, SF-12, and VAS scores, as well as reduced rates of secondary surgeries and clinically relevant adjacent segment disease (ASD), defined as at least one increase in Kellgren-Lawrence grade from baseline. At 2 years postop, radiographs showed Mobi-C to support improved mean ROM both flexion-extension ( $10.8^{\circ} \pm 6.5^{\circ}$ ) and lateral bending ( $5.4^{\circ} \pm 3.1^{\circ}$ ). Data collection from the Mobi-C IDE study has continued since the initial 2-year data was published, with 7-year data continuing to support the safety and efficacy of Mobi-C.<sup>5</sup>

**Table 3. 2-Year Results of 1-Level Mobi-C<sup>1,13</sup>**

	Mobi-C	ACDF
<b>Follow-up</b>	94.3%	92.0%
<b>Success</b>	No serious device- or device procedure-related AEs; no supplemental surgical procedure at the index level; maintenance or improvement of neurological function compared to baseline; and $\geq 15$ point improvement in NDI (on 100-point scale)	
<b>Definition</b>		
<b>Success Rate</b>	73.7%	65.3%
<b>Mobi-C is non-inferior to ACDF at 1-level</b>		
<b>Progression of ASD</b>		
<b>Superior</b>	15.0%	26.1%
<b>Inferior</b>	8.3%	21.0%
<b>Clinically Relevant ASD</b>		
<b>Superior</b>	9.1%	14.4%
<b>Inferior</b>	2.7%	11.1%
<b>Heterotopic Ossification</b>		
<b>Grade 3</b>	9.5%	-
<b>Grade 4</b>	5.9%	-
<b>Grade 3/4 Combined</b>	15.4%	-
<b>Secondary Surgery</b>		
<b>Index level</b>	1.2%	6.2%
<b>Adjacent level</b>	1.1%	3.7%
<b>NDI Score</b>		
<b>Pre-op</b>	53.4	54.1
<b>2-Year</b>	16.3	19.9
<b>NDI Improvement</b>	37.1	34.2
<b>Range of Motion (ROM)</b>		
<b>Flexion-Extension</b>	10.7	-
<b>Lateral Bending</b>	5.4	-

While direct clinical comparison between Mobi-C and M6-C can't be made due to the lack of a head-to-head clinical study, the study design and data reporting in the M6-C Summary of Safety and Effectiveness (SSED)<sup>2</sup> leave something to be desired. The study was non-randomized, with site-specific treatment groups, meaning that no surgeons in the study performed both cervical total disc replacement (TDR) and ACDF. The ACDF control in the study was comprised of two cohorts – one was treated concurrently with the TDR treatment group at different centers, while the other consisted of historical ACDF control patients from a previous cervical disc study. Due to differences in data collection and outcome measures between the concurrent and historical patients, some data from the control cohorts could not be pooled via a validated method. This prevented direct comparison of certain outcomes (SF-12/SF-36, Odom's Criteria, and Patient Satisfaction) between M6-C and ACDF. When direct comparisons were possible, the M6-C cohort had higher success rates, greater mean NDI score improvement, and lower rates of secondary surgery than the ACDF cohort.

**Table 4. 2-Year Results of 1-Level M6-C<sup>2</sup>**

	M6-C	ACDF
<b>Follow-up</b>	95.0%	87.7%
<b>Success</b>	No study failure due to secondary surgical interventions at the index level; absence of major complications defined as radiographic failure, neurologic failure, or failure by AE; and NDI score improvement of $\geq 15/50$ points in subjects with a baseline score $\geq 30/50$ points, or a 50% improvement in subjects with a baseline score $< 30/50$ points	
<b>Definition</b>		
<b>Success Rate</b>	86.8%	79.3%
<b>M6-C is non-inferior to ACDF at 1-level</b>		
<b>Progression of ASD</b>		
<b>Superior</b>	NR	NR
<b>Inferior</b>	NR	NR
<b>Clinically Relevant ASD</b>		
<b>Superior</b>	NR	NR
<b>Inferior</b>	NR	NR
<b>Heterotopic Ossification</b>		
<b>Grade 3</b>	10.7%	-
<b>Grade 4</b>	0.7%	-
<b>Grade 3/4 Combined</b>	11.3%	-
<b>Secondary Surgery</b>		
<b>Index level</b>	1.9%	4.8%
<b>Adjacent level</b>	3.1%	2.1%
<b>NDI Score</b>		
<b>Pre-op</b>	54.8	51.9
<b>2-Year</b>	12.1	17.9
<b>NDI Improvement</b>	42.7	34.0
<b>Range of Motion (ROM)</b>		
<b>Flexion-Extension</b>	8.8	-
<b>Lateral Bending</b>	6.9	-

## References:

1. Hisey MS, Bae HW, Davis R, et al. Multi-center, Prospective, Randomized, Controlled Investigational Device Exemption Clinical Trial Comparing Mobi-C Cervical Artificial Disc to Anterior Discectomy and Fusion in the Treatment of Symptomatic Degenerative Disc Disease in the Cervical Spine. *Int J Spine Surg* 2014;8:7.
2. U.S. Food and Drug Administration. "M6-C® Summary of Safety and Effectiveness Data." [https://www.accessdata.fda.gov/cdrh\\_docs/pdf17/P170036B.pdf](https://www.accessdata.fda.gov/cdrh_docs/pdf17/P170036B.pdf). Accessed on May 1, 2019.
3. Jadik S, Miller FC, Pietila AT. Cervical arthroplasty with the M6-C artificial disc in degenerative disc disease of the cervical spine. *Arch Med* 2015;7(5).
4. Laurysen C, Coric D, Dimmig T, et al. Cervical total disc replacement using a novel compressible prosthesis: Results from a prospective Food and Drug Administration–regulated feasibility study with 24-month follow-up. *Int J Spine Surg* 2012;6:71-77.
5. Radcliff K, Davis RJ, Hisey MS, et al. Long-term evaluation of cervical disc arthroplasty with the Mobi-C cervical disc: A randomized, prospective, multicenter trial with seven-year follow-up. *Int J Spine Surg* 2017;11(31).
6. White AA III, Panjabi MM. Kinematics of the spine. In: *Clinical Biomechanics of the Spine*, 2nd ed. Baltimore, MD: Lippincott Williams & Wilkins; 1990:97-98.
7. Liu J, Ebraheim NA, Haman SP, et al. How the Increase of the Cervical Disc Space Height Affects the Facet Joint: An Anatomy Study. *Spine* 2006;31(12):E350-E354.
8. Peng CWB, Quirnoa M, Bendo JA, et al. Effect of intervertebral disc height on postoperative motion and clinical outcomes after Prodisc-C cervical disc replacement. *Spine J* 2009;9(7):551-555.
9. Yuan W, Zhang H, Zhou X, et al. The Influence of Artificial Cervical Disc Prosthesis Height on the Cervical Biomechanics: A Finite Element Study. *World Neurosurg* 2018;113:e490-e498.
10. Staudt MD, Das K, Duggal N. Does design matter? Cervical disc replacements under review. *Neurosurg Rev* 2018;41(2):399-407.
11. Kim KD, Karnati T, Permeswaran VJ, et al. Overstuffing a total disc replacement device significantly reduces index level and total range of motion: a biomechanical study [Internal white paper].
12. Park CK, Ryu KS. Are Controversial Issues in Cervical Total Disc Replacement Resolved or Unresolved?: A Review of Literature and Recent Updates. *Asian Spine J* 2018;12(1):178-192.
13. Mobi-C IDE Study 2-year data on file, available upon request.

