

Quality and Confidence –  
All in One

# Screw Fixation System



## 1. Secure and Simple

- Easy and stable screw pick-up
- Reliable driver to screw connection
- Stored and sterilized in one durable kit tray

## 2. Compact and Complete

- Autoclavable metal storage tray
- Heavy, stainless steel screw driver handle
- 1.5 mmD Micro and 2.0 mmD Mini System

## 3. Versatile Module – Mix and Match

- Titanium alloy fixation screws in 1.5 mmD, 2.0 mmD and 2.3 mmD
- Color coded drills, drivers and screws
- Choose and add for your individual practice needs at any time

 **ZimVie**

# Quality and Confidence – All in One

The Screw Fixation System provides a compact solution for the temporary fixation and stabilization of bone transplants, suitable resorbable and non-resorbable bone replacement materials, and membranes for the alveolar ridge. Two color coded systems in 1.5 mmD Micro (Blue) and 2.0 mmD Mini (Red) offer concise and cost effective functional options. The color coding scheme for the two systems, the components and the screws makes rapid identification of the parts possible and simplifies parts matching. This modular storage system permits individual configuration and the open design ensures optimum access during cleaning and sterilization. Fixation screws and mesh are manufactured from pure titanium or titanium alloy. They are biocompatible, corrosion-proof and nontoxic in the biological environment. They allow imaging virtually free of artifacts.



## Secure and Simple

The screws are easily picked up and have a reliable connection to the driver.

## Getting Started Made Easy with Our Start-Up Screw Fixation Kit

### Start-Up Screw Fixation Kit\* (69.01.10Z ) features:

- Durable, autoclavable metal storage tray
- Heavy stainless steel screw driver handle
- Micro System, 1.5 mmD
  - Stainless steel screw driver inserts (long and short)
  - Stainless steel pilot drill, 14 mmL
  - Stainless steel pilot block drill, 15 mmL



## Screw Fixation kit – Start-Up Kit Ordering Information

Contents	Tray	Screw Driver Handle	Screw Driver Insert, Short	Screw Driver Insert, Long	Pilot Drill, Micro, 14 mmL	Pilot Block Drill, Micro
Assembled Start-Up Kit*	69.01.10Z					
Components	69.01.11Z	75.23.52Z	75.23.23Z	75.23.19Z	69.01.09Z	69.01.16Z

\*Screws are not included in the start-up kit

For more information, visit [ZimVie.com](http://ZimVie.com)

**ZimVie**  
 4555 Riverside Drive  
 Palm Beach Gardens, FL 33410  
 1-800-342-5454

**ZimVie**  
 Restoring Daily Life.®