

An Optimal Blend of Scaffold and Osteoinductive Potential



Puros® Ci
Particulate Allograft



Ideal Blend¹

- Combination of mineralized and demineralized allograft bone (DBM)
- Only DBM with verified osteoinductive (OI) potential* is released for processing
- Mineralized component acts as a scaffold; Mineralized particles include both cancellous and cortical bone chips

Best In Class Processing¹

- Proprietary Tutoplast® processing preserves the native collagen matrix of the mineralized particles
- Proprietary Cancelli SP® processing preserves OI potential

Trusted Source¹

- Safety and reliability you expect from the Puros Brand
- Single donor source for mineralized and demineralized particles
- Terminally sterilized to SAL 10⁻⁶

Convenient Handling

- Quick hydration
- Room temperature storage²



*DBM is evaluated for bone formation utilizing an in-vivo athymic rat model. Findings from an animal model are not necessarily predictive of human clinical results.

Puros® Ci Particulate Allograft

A unique combination of cancellous, cortical and demineralized cortical bone designed for filling bone voids.


Tutoplast Process¹

The Tutoplast Process is a proprietary process that preserves biomechanical/biochemical integrity and native tissue structure. Over the past 24 years:

- Approximately 7 million implants have gone through the Tutoplast Process globally
- More than 5.5 million patient procedures have been performed using implants processed through the Tutoplast Process globally
- Zero confirmed incidence of implant-associated infection

Cancellle SP process¹

The Cancellle SP Process is a proprietary process that sterilizes DBM while simultaneously allowing it to maintain its OI potential.



Ordering Information:

Item Number	Description
69800	Puros Ci Particulate Allograft, 0.5 cc, .25 – 1 mm
69801	Puros Ci Particulate Allograft, 1 cc, .25 – 1 mm
69802	Puros Ci Particulate Allograft, 2 cc, .25 – 1 mm
69803	Puros Ci Particulate Allograft, 0.5 cc, 1 – 2 mm
69804	Puros Ci Particulate Allograft, 1 cc, 1 – 2 mm
69805	Puros Ci Particulate Allograft, 2 cc, 1 – 2 mm

References

1. Data on file with RTI Surgical, Inc.
2. Refer to product label for specific storage parameters.

For more information, visit ZimVie.com

ZimVie
4555 Riverside Drive
Palm Beach Gardens, FL 33410
1-800-342-5454
Phone: +1-561-776-6700
Fax: +1-561-776-1272

 **ZimVie**
Restoring Daily Life.®

Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. Cancellle SP is a U.S. registered trademark of RTI Surgical, Inc. Puros is manufactured by RTI Surgical, Inc. Tutoplast is a registered trademark of Tutogen Medical GmbH. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV1467 REV A 07/23 ©2023 ZimVie Inc. All rights reserved.

