




Kit Assembly

T3[®] Short Implant System Surgical Kit



 **ZimVie**

ZimVie DENTAL SOLUTIONS

Kit Assembly - T3 Short Implant System Surgical Kit (BSKIT)



1

ACT®
Pointed Starter
Drill ACTPSD



2

ACT Drill
2 mmD x
5-6 mmL
ACT206S




3

ACT Drill
3.25 mmD x
5-6 mmL
ACT326S



4

ACT Drill
3.85 mmD x
5-6 mmL
ACT386S



4

ACT Drill
4.25 mmD x
5-6 mmL
ACT426S



5

ACT Drill
4.85 mmD x
5-6 mmL
ACT486S



5S

Shaping
Drill
3.85 mmD x
5 mmL
FCS385S



5L

Shaping Drill
3.85 mmD x
6 mmL
FCS386S



6S

Shaping Drill
4.85 mmD x
5 mmL
FCS485S



6L

Shaping Drill
4.85 mmD x
6 mmL
FCS486S



7

Dense Bone
Tap
5 mmD x
5-6 mmL
TAP56S



7

Dense Bone
Tap
6 mmD x
5-6 mmL
TAP66S



8

Handpiece
Connector
MDR10



9

Large Hex
Driver
PHD02N




10

Short Ratchet
Extension
RE100



11

Small Hex
Driver
PHD00N



12

High Torque
Indicating Ratchet
Wrench
H-TIRW



13

Open End
Wrench
CW100

For more information, visit ZimVie.com

ZimVie
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272
dentalCS@ZimVie.com

ZimVie
Restoring Daily Life.®

Unless otherwise indicated, as referenced herein, all trademarks are the property of ZimVie; and all products are manufactured by one or more of the dental subsidiaries of ZimVie, Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV1562 REV A 09/23 ©2023 ZimVie. All rights reserved.

