



Touchless transfer.  
Superior control.



# Internal Connection Universal Placement Driver Tips

With GemLock® Technology for  
the Certain® Connection



1

## Secure Fit

Touchless delivery tools offer a secure, reliable fit

2

## Consistent Performance

Durable synthetic rubies improve driver tip longevity and ensure delivery tools remain consistent over multiple uses

3

## Superior Control

GemLock Technology engages components a for precise handling



# Internal Connection Universal Placement Driver Tips

With GemLock Technology for the Certain Connection

## Secure Fit

- Internally engaging driver tips that attach to surgical handpieces in order to pick up and carry implants, healing abutments, and cover screws with the Certain Connection
- No drop delivery from package to patient increases efficiency, decreases chair time, and minimizes the potential for contamination of the implant

## Consistent Performance

- Spring-loaded synthetic ruby bearings seamlessly and repeatedly secure the implant to the delivery tool
- Synthetic rubies have extreme abrasion resistance and low breakaway friction resulting in expanded longevity and performance of your driver tips

## Superior Control

- Superior delivery with dynamic tactile response and touchless transfer during implant surgery for a simpler, more controlled experience
- Engages components internally for increased precision, making them ideal for use in tight interdental spaces



GemLock Technology

## Ordering Information

Item Number	Description
IIPDTUL	Internal Connection Universal Placement Driver Tip (Long)
IIPDTUS	Internal Connection Universal Placement Driver Tip (Short)

For more information, visit [ZimVie.com](https://www.zimvie.com)

**ZimVie**  
4555 Riverside Drive  
Palm Beach Gardens, FL 33410  
1-800-342-5454

 **ZimVie**  
Restoring Daily Life.®

Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV1912 REV A 12/23 ©2023 ZimVie. All rights reserved.

