

# RealGUIDE™ Z3D Guided Surgery Kit

Designed for Zimmer Biomet's TSV®, Trabecular Metal™ and Eztetic® Dental Implants

FULLY  
GUIDED!



 RealGUIDE™ Z3D  
GUIDED SURGERY KIT



## RealGUIDE Z3D GUIDED SURGERY KIT

The new RealGUIDE Z3D Guided Surgery Kit contains instruments needed for the planned surgery. It is designed to work with Zimmer Biomet's **TSV, Trabecular Metal and Eztetic** Implant Systems.

The kit information is available in the RealGUIDE Software Suite, an open system for diagnosis and advanced surgical guide design that will support you with customized drilling reports to be well-prepared for dental implant therapy.

- Fully guided placement of Zimmer Biomet Dental Eztetic 3.1 mm Implant and TSV/TM Implants in sizes 3.7 mm, 4.1 mm and 4.7 mm
- Color-coded for ease of use
- Compact kit, which is useful for space-constrained dental offices
- No handles or keys are necessary while preparing the osteotomy
- DLC (Diamond-like Carbon) coating on the drills provides wear resistance characteristics
- Includes implant mounts hex feature to align with the hex feature on the corresponding guide sleeves
- Staging pedestal for holding the implant vial to remove the Fixture Mount Transfer (FMT) of the TSV implants is available. This enables the TSV implants to be placed fully guided through the surgical guides
- Both soft and dense bone protocols are available



RealGUIDE Software Suite  
with PLAN, APP, GUIDE Modules

Product	Part No.
RealGUIDE Z3D Surgical Kit	3DM0070ZUS



MORE INFORMATION  
ABOUT THE RealGUIDE Z3D  
GUIDED SURGERY KIT AT  
[ZBDGUIDEDSURGERY.COM](http://ZBDGUIDEDSURGERY.COM)

### Contact us

at 1-800-342-5454 or visit [zbdguidedurgery.com](http://zbdguidedurgery.com)



Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc., and distributed and marketed by Zimmer Biomet Dental (and, in the case of distribution and marketing, its authorized marketing partners). The RealGUIDE Z3D is manufactured by 3Diemme, a subsidiary of Zimmer Biomet Holdings. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB1300 REV A 04/21 ©2021 Zimmer Biomet. All rights reserved.

