

RealGUIDE™ CAD Software

All-inclusive, flexible, and powerful dental design



RealGUIDE™
UNIVERSAL OPEN SYSTEM



ALL-INCLUSIVE DENTAL DESIGN SOFTWARE


The RealGUIDE™ CAD Software is a powerful, all-in-one dental software for restorative design. As part of the RealGUIDE Full Suite it offers virtually boundless possibilities – from simple to highly complex dental designs without the need to purchase any separate modules. The all-in-one software solution offers you maximum flexibility and all functionalities you might need at any time – with no hidden costs.



Get a free trial of the RealGUIDE™ Full Suite now for 30 days! Get more information at realguide.com

FEATURES

- Free-form tools and trademarked SandBox™ feature for complex restorations
- Flows organized in wizards with a cascade structure and icons
- Single-step modelling of removable dentures and bars
- One single software immediate loading protocol
- Share files with your partners on the multidisciplinary team



Unprecedented Freedom and flexibility

Flows organized in wizards with a cascade structure and icons that allow you to start with a generic reconstruction type and specify the details during design. Free-form tools give you the flexibility and freedom to work in the same way as traditional dentistry. Simplified workflows support single-step modelling of removable dentures and bars based on specific libraries. And with the cloud licensing instead of dongles, realize efficiencies in your day-to-day operations.


INDICATIONS

- Inlay, Onlay, Richmond, Veneers
- Crowns & Bridges
- Individual and telescopic abutments
- Bars and telescopic bars
- Overdenture and removable prosthesis
- Full Dentures
- Immediate Load Provisional
- Pre-OP
- Models and orthodontic Models
- Bite Splints
- Digital Waxups
- Thimbles
- Attachments

FUNCTIONALITIES

- Provisional Planning
- Virtual Articulator
- Individual TMJ Movements
- Jaw Motion Visualizer
- Virtual parallelometer
- Bar Designer
- Model Creator
- Bite Splint Module

FREE
FOR RealGUIDE
SUBSCRIBERS



SandBox™ FEATURE

The trademarked SandBox™ is a set of tools and processes that allow you to create complex cases and build customized workflows. Ensure data integrity and time savings without the need to go to another software.

- Create customized CAD environment with individual protocols and libraries
- Generate any geometric primitives – no additional software needed
- Cut, split, merge, repair any STL file in a single, integrated environment
- Patented Volume-Splitting-Function

SOFTWARE FEATURE OVERVIEW

	 PLAN	 APP	 GUIDE	 CAD	 FULL SUITE
	FREE!		 INCLUDED		30-DAY FREE TRIAL
Free implant treatment planning for life	✓	×	✓	×	✓
Vast implant library	✓	×	✓	✓	✓
Limited Artificial Intelligence tools	✓	×	✓	✓	✓
Extensive Artificial Intelligence tools	×	×	✓	✓	✓
Guide design and print	×	×	✓	×	✓
Affordable monthly subscription and click fees	×	✓	✓	✓	✓
Export STL files optimized for 3D printing or milling	×	×	✓	✓	✓
Message and share data with other team members	✓	✓	✓	✓	✓
Download the RealGUIDE™ App to your iOS devices	×	✓	✓	✓	✓
Manage your projects from any location**	✓	✓	✓	✓	✓
View and superimpose files, trace anatomical structures and place implants with easy clicks	✓	×	✓	×	✓
View plans and tweak implant positions with your fingertip or with easy clicks	✓	✓	✓	✓	✓
Restorative Design and manufacturing from simple to complex constructions	×	×	×	✓	✓
Immediate load provisional prosthesis modelling	×	×	×	✓	✓
Provisional Planning, Virtual Articulator, TMJ Movements, Jaw Motion Visualizer, Virtual Parallelometer, Bar Designer, Model Creator, Bite Splint module	×	×	×	✓	✓
Free-form tools and patented Sandbox Feature (for unusual cases and stackable guides)	×	×	✓	✓	✓
Max. number of patients in the cloud	5	5/100 ***	100	100	100

* Available as an option. ** As long as the device is connected on the internet and has the appropriate software installed.

*** Depending on the installed software module "PLAN" or "GUIDE".

Reference:

Virtual Implant Rehabilitation of the Severely Atrophic Maxilla: A Radiographic Study, Michele Manacorda, Bianca Poletti de Chaurand, Alberto Merlone, Giulia Tetè, Francesca Mottola and Raffaele Vinci Dent J (Basel). 2020 Mar; 8(1): 14. Published online 2020 Feb 2. doi: 10.3390/dj8010014



MORE INFORMATION ABOUT THE DIFFERENT MODULES AT [REALGUIDE.COM](https://realguide.com)



Contact us at 1-866-828-5300 or visit realguide.com

