



The Immediacy
Patients Want





Certain[®] Tapered Navigator[®] System For Guided Surgery

Kit Assembly



Kit Assembly Layout






- 
40
SGTTP4
- 
41
SGTTP5
- 
42
SGTCP4S
- 
43
SGTCP5S
- 
44
SGTCP4L
- 
45
SGTCP5L
- 
46
SGTD19A
- 
47
SGTD19B
- 
48
SGTD19C
- 
49
SGTD19D
- 
50
SGTD19E
- 
51
SGTD19F
- 
52
SGTD19G
- 
53
MSGTCD5
- 
54
SGTCD4S
- 
55
SGTCD5S
- 
56
MSGTCDL
- 
57
SGTCD4L
- 
58
SGTCD5L
- 
59
MSGD85S
- 
60
MSGTD10S
- 
61
MSGTD11S
- 
62
MSGTD13S
- 
63
MSGTD15S
- 
63
SGTD485S
- 
64
SGTD410S
- 
65
SGTD411S
- 
78
MSGTD15L
- 
79
SGTD485L
- 
80
SGTD410L
- 
81
SGTD411L
- 
82
SGTD413L
- 
83
SGTD415L
- 
84
SGTD585L
- 
85
SGTD510L
- 
86
SGTD511L
- 
87
SGTD513L
- 
88
SGTD515L



Implant Mount Assembly Layout




39


SGTTRAY



Tap Kit Assembly Layout



Implant Analog Mount Kit



1
SGLTRAY

Tapered Navigator System Instrumentation

	Item No.	Description
1	SGTLTRAY	Tapered Navigator Laboratory Tray
2	MSGIAMS	Navigator Certain® Analog Mount, 3.4 mm(D) - Short
3	MSGIAML	Navigator Certain Analog Mount, 3.4 mm(D) - Long
4	SGIAM4S	Navigator Certain Analog Mount, 4.1 mm(D) - Short
5	SGIAM4L	Navigator Certain Analog Mount, 4.1 mm(D) - Long
6	SGIAM5S	Navigator Certain Analog Mount, 5 mm(D) - Short
7	SGIAM5L	Navigator Certain Analog Mount, 5 mm(D) - Long
8	SGTTAPT	Tapered Navigator Bone Tap Tray
9	MSGTT85S	Tapered Navigator Bone Tap, 3.25 mm(D) x 8.5 mm(L) - Short
10	MSGTT10S	Tapered Navigator Bone Tap, 3.25 mm(D) x 10 mm(L) - Short
11	MSGTT11S	Tapered Navigator Bone Tap, 3.25 mm(D) x 11.5 mm(L) - Short
12	MSGTT13S	Tapered Navigator Bone Tap, 3.25 mm(D) x 13 mm(L) - Short
13	MSGTT15S	Tapered Navigator Bone Tap, 3.25 mm(D) x 15 mm(L) - Short
14	SGTT485S	Tapered Navigator Bone Tap, 4 mm(D) x 8.5 mm(L) - Short
15	SGTT410S	Tapered Navigator Bone Tap, 4 mm(D) x 10 mm(L) - Short
16	SGTT411S	Tapered Navigator Bone Tap, 4 mm(D) x 11.5 mm(L) - Short
17	SGTT413S	Tapered Navigator Bone Tap, 4 mm(D) x 13 mm(L) - Short
18	SGTT415S	Tapered Navigator Bone Tap, 4 mm(D) x 15 mm(L) - Short
19	SGTT585S	Tapered Navigator Bone Tap, 5 mm(D) x 8.5 mm(L) - Short
20	SGTT510S	Tapered Navigator Bone Tap, 5 mm(D) x 10 mm(L) - Short
21	SGTT511S	Tapered Navigator Bone Tap, 5 mm(D) x 11.5 mm(L) - Short
22	SGTT513S	Tapered Navigator Bone Tap, 5 mm(D) x 13 mm(L) - Short
23	SGTT515S	Tapered Navigator Bone Tap, 5 mm(D) x 15 mm(L) - Short
24	MSGTT85L	Tapered Navigator Bone Tap, 3.25 mm(D) x 8.5 mm(L) - Long
25	MSGTT10L	Tapered Navigator Bone Tap, 3.25 mm(D) x 10 mm(L) - Long
26	MSGTT11L	Tapered Navigator Bone Tap, 3.25 mm(D) x 11.5 mm(L) - Long
27	MSGTT13L	Tapered Navigator Bone Tap, 3.25 mm(D) x 13 mm(L) - Long
28	MSGTT15L	Tapered Navigator Bone Tap, 3.25 mm(D) x 15 mm(L) - Long
29	SGTT485L	Tapered Navigator Bone Tap, 4 mm(D) x 8.5 mm(L) - Long
30	SGTT410L	Tapered Navigator Bone Tap, 4 mm(D) x 10 mm(L) - Long
31	SGTT411L	Tapered Navigator Bone Tap, 4 mm(D) x 11.5 mm(L) - Long
32	SGTT413L	Tapered Navigator Bone Tap, 4 mm(D) x 13 mm(L) - Long
33	SGTT415L	Tapered Navigator Bone Tap, 4 mm(D) x 15 mm(L) - Long
34	SGTT585L	Tapered Navigator Bone Tap, 5 mm(D) x 8.5 mm(L) - Long
35	SGTT510L	Tapered Navigator Bone Tap, 5 mm(D) x 10 mm(L) - Long
36	SGTT511L	Tapered Navigator Bone Tap, 5 mm(D) x 11.5 mm(L) - Long
37	SGTT513L	Tapered Navigator Bone Tap, 5 mm(D) x 13 mm(L) - Long
38	SGTT515L	Tapered Navigator Bone Tap, 5 mm(D) x 15 mm(L) - Long
39	SGTTRAY	Tapered Navigator Surgical Tray
40	SGTTP4	Tapered Navigator Tissue Punch, 4.1 mm(D)
41	SGTTP5	Tapered Navigator Tissue Punch, 5 mm(D)
42	SGTCP4S	Tapered Navigator Cortical Perforator, 4.1 mm(D) - Short
43	SGTCP5S	Tapered Navigator Cortical Perforator, 5 mm(D) - Short
44	SGTCP4L	Tapered Navigator Cortical Perforator, 4.1 mm(D) - Long
45	SGTCP5L	Tapered Navigator Cortical Perforator, 5 mm(D) - Long
46	SGTD19A	Tapered Navigator Twist Drill, 1.9 mm(D) x a(L)
47	SGTD19B	Tapered Navigator Twist Drill, 1.9 mm(D) x b(L)
48	SGTD19C	Tapered Navigator Twist Drill, 1.9 mm(D) x c(L)
49	SGTD19D	Tapered Navigator Twist Drill, 1.9 mm(D) x d(L)
50	SGTD19E	Tapered Navigator Twist Drill, 1.9 mm(D) x e(L)
51	SGTD19F	Tapered Navigator Twist Drill, 1.9 mm(D) x f(L)
52	SGTD19G	Tapered Navigator Twist Drill, 1.9 mm(D) x g(L)
53	MSGTCDS	Tapered Navigator Countersink Drill, 3.4 mm(D) - Short

	Item No.	Description
54	SGTCD4S	Tapered Navigator Countersink Drill, 4.1 mm(D) - Short
55	SGTCD5S	Tapered Navigator Countersink Drill, 5 mm(D) - Short
56	MSGTCDL	Tapered Navigator Countersink Drill, 3.4 mm(D) - Long
57	SGTCD4L	Tapered Navigator Countersink Drill, 4.1 mm(D) - Long
58	SGTCD5L	Tapered Navigator Countersink Drill, 5 mm(D) - Long
59	MSGTD85S	Tapered Navigator Shaping Drill, 3.25 mm(D) x 8.5 mm(L) - Short
60	MSGTD10S	Tapered Navigator Shaping Drill, 3.25 mm(D) x 10 mm(L) - Short
61	MSGTD11S	Tapered Navigator Shaping Drill, 3.25 mm(D) x 11.5 mm(L) - Short
62	MSGTD13S	Tapered Navigator Shaping Drill, 3.25 mm(D) x 13 mm(L) - Short
63	MSGTD15S	Tapered Navigator Shaping Drill, 3.25 mm(D) x 15 mm(L) - Short
64	SGTD485S	Tapered Navigator Shaping Drill, 4 mm(D) x 8.5 mm(L) - Short
65	SGTD410S	Tapered Navigator Shaping Drill, 4 mm(D) x 10 mm(L) - Short
66	SGTD411S	Tapered Navigator Shaping Drill, 4 mm(D) x 11.5 mm(L) - Short
67	SGTD413S	Tapered Navigator Shaping Drill, 4 mm(D) x 13 mm(L) - Short
68	SGTD415S	Tapered Navigator Shaping Drill, 4 mm(D) x 15 mm(L) - Short
69	SGTD585S	Tapered Navigator Shaping Drill, 5 mm(D) x 8.5 mm(L) - Short
70	SGTD510S	Tapered Navigator Shaping Drill, 5 mm(D) x 10 mm(L) - Short
71	SGTD511S	Tapered Navigator Shaping Drill, 5 mm(D) x 11.5 mm(L) - Short
72	SGTD513S	Tapered Navigator Shaping Drill, 5 mm(D) x 13 mm(L) - Short
73	SGTD515S	Tapered Navigator Shaping Drill, 5 mm(D) x 15 mm(L) - Short
74	MSGTD85L	Tapered Navigator Shaping Drill, 3.25 mm(D) x 8.5 mm(L) - Long
75	MSGTD10L	Tapered Navigator Shaping Drill, 3.25 mm(D) x 10 mm(L) - Long
76	MSGTD11L	Tapered Navigator Shaping Drill, 3.25 mm(D) x 11.5 mm(L) - Long
77	MSGTD13L	Tapered Navigator Shaping Drill, 3.25 mm(D) x 13 mm(L) - Long
78	MSGTD15L	Tapered Navigator Shaping Drill, 3.25 mm(D) x 15 mm(L) - Long
79	SGTD485L	Tapered Navigator Shaping Drill, 4 mm(D) x 8.5 mm(L) - Long
80	SGTD410L	Tapered Navigator Shaping Drill, 4 mm(D) x 10 mm(L) - Long
81	SGTD411L	Tapered Navigator Shaping Drill, 4 mm(D) x 11.5 mm(L) - Long
82	SGTD413L	Tapered Navigator Shaping Drill, 4 mm(D) x 13 mm(L) - Long
83	SGTD415L	Tapered Navigator Shaping Drill, 4 mm(D) x 15 mm(L) - Long
84	SGTD585L	Tapered Navigator Shaping Drill, 5 mm(D) x 8.5 mm(L) - Long
85	SGTD510L	Tapered Navigator Shaping Drill, 5 mm(D) x 10 mm(L) - Long
86	SGTD511L	Tapered Navigator Shaping Drill, 5 mm(D) x 11.5 mm(L) - Long
87	SGTD513L	Tapered Navigator Shaping Drill, 5 mm(D) x 13 mm(L) - Long
88	SGTD515L	Tapered Navigator Shaping Drill, 5 mm(D) x 15 mm(L) - Long
89	MDR10	Handpiece Connector
90	PHD02N	Narrow Posterior Large Hex Driver, 17 mm
91	RASH3N	Narrow Right Angle Large Hexed Driver Tip, 24 mm
92	RE100	Ratchet Extension 6 mm(L)
93	RE200	Ratchet Extension 15 mm(L)
94	SGTDPHT	Tapered Navigator Twist Drill Positioning Handle
95	SGTDPHR	Tapered Navigator Reduction Drill Positioning Handle
96	CW100	Open End Wrench
97	WR150	Ratchet Wrench
98	MSGIIMS	Navigator Certain Implant Mount, 3.4 mm(D) - Short
99	MSGIIML	Navigator Certain Implant Mount, 3.4 mm(D) - Long
100	SGIIM4S	Navigator Certain Implant Mount, 4.1 mm(D) - Short
101	SGIIM4L	Navigator Certain Implant Mount, 4.1 mm(D) - Long
102	SGIIM5S	Navigator Certain Implant Mount, 5 mm(D) - Short
103	SGIIM5L	Navigator Certain Implant Mount, 5 mm(D) - Long
104	IMMBP	Certain Manual Bone Profiler, 3.4 mm(D) x 4.1 mm(P)
105	IMBP4	Certain Manual Bone Profiler, 4.1 mm(D) x 5 mm(P)
106	IMBP5	Certain Manual Bone Profiler, 5 mm(D) x 6 mm(P)

For more information, visit [ZimVie.com](https://www.zimvie.com)

ZimVie
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272
dentalCS@ZimVie.com



Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV1210 REV A 05/23 ©2023 ZimVie. All rights reserved.

