

Your path to predictable

Unleash an exceptional digital workflow and scanning experience that drives clinic efficiency, restorative solutions, and patient and dentist comfort

Powerful hardware and intelligent software, the full-feature Medit i700 is the key to unlocking your clinic's full potential.

Medit i700 scanner

Plug and Scan Directly plugs into a PC with a single cable
Wireless option available!





Medit Link Secure, cloud-based 3D data sharing platform that facilitates communication and collaboration with other dental partners
Compatible with numerous systems, including **RealGUIDE®***, for enhanced workflows.
*anticipated 2024

Medit Apps A wide range of features to support restorative outcomes:

- Scan and color filtering
- Margin line creation
- Ortho simulation
- Temporaries and crown fit
- A.I. scan body matching



Unleash your clinic's full potential

Category		MEDIT i600	MEDIT i700	MEDIT i700 wireless	MEDIT i900
Model					
Scanning Technology	Capture Rate	Fast	Faster	Faster	Fastest
	Imaging Technology	3D-in-motion video technology/3D full color streaming capture			
	Anti-Fogging Technology	Adaptive air flow			
Accuracy	Full Arch	10.9 um + 0.98			
Handpiece	Dimensions	248 x 44 x 47.4 mm		313 x 44 x 47.4 mm	223 X 37 X 35 mm
	Total Weight	245 g		328 g	165 g
Tip	Scan Area	15 x 13 mm			18 x 15 mm
	Autoclavable	150 times			Sterilize for 10 minutes at 135°C (275°F) at gravity type and dry for 30 minutes. Sterilize for 4 minutes at 134°C (273°F) at pre-vacuum type and dry for 20 minutes.
		Autoclave 121° C 30 mins (Gravity Type) Autoclave 134° C 4 mins (Pre-Vacuum Type) Autoclave 135° C 10 min (Gravity Type)			
	Reversable Tip	Yes	Yes	Yes	No
Special Mode	Remote Control Mode	No	Yes		
	UV-LED Disinfection	No	Yes		
	Smart Shade Guide	Yes			
	HD Mode	No	Yes		
MAC OS Compatibility		Yes			
Connectivity	Connection Type	Wired USB 3.1 Gen 1		Wireless USB 3.1 Gen 1	Wired USB 3.0
	Plug and Scan	Yes			
	Detachable Cable	Yes			

For more information, visit ZimVie.com

ZimVie
4555 Riverside Drive
Palm Beach Gardens, FL 33410
1-800-342-5454
Phone: +1-561-776-6700
Fax: +1-561-776-1272



ZimVie Education



Medit User Group



Medit Academy



Medit Help Center

 **ZimVie**
Restoring Daily Life.®

Unless otherwise indicated, as referenced herein, all trademarks are the property of ZimVie; or and affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. Intraoral scanners and accessories are manufactured by Medit Corp. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV2430 GLBL REV A 06/24 ©2024 ZimVie. All rights reserved.

