

Restorative Tooling Chart

Certain®

Torque 20 Ncm | .048/1.20 mm Hex

Option 1

1



L-TIRW

IHX1.2D/IXL1.2D* or RASH3N/RASH8N

Option 2

2



TWR with ARA-DTSQ (Azure™ Adapter)

IHX1.2/IXL1.2D* or RASH3N/RASH8N


Restorative Driver and Tips

Length	Finger Driver	Driver Tip	Finger Driver	Driver Tip	Finger Driver	Driver Tip	Driver Tip	Finger Driver	Driver Tip
17 mm	PAD00		PHD02N	IHX1.2D*/ RASH3N	ICD00				
20 mm								HXGR1.25*	
21 mm							ZFX02HLD21		
23 mm									HX1.25D*
24 mm	PAD24	RASA3	PHD03N	IXL1.2D*/ RASH8N		ASCDT24	TASCT24		
26 mm									HXL1.25D*
28 mm							ZFX02HLD28		
30 mm						ASCDT30	TASCT30	HXLGR1.25*	
36 mm							ZFX02HLD36		
For use with:	Straight Low Profile Abutments		Healing Abutments, Abutments, Prosthetic Screws, Cover Screws (.048/1.20 mm hex) Angled Low Profile Abutments		Twist Lock Transfer Impression Coping	BellaTek® Flex and Express Ti-Bases	GenTek® Angulated Screw Channel Ti-Bases	GenTek Scan Bodies (.050/1.25 mm hex)	

*GemLock®

Torque Wrenches and Adapters

Torque Wrench	Item No.
Low Torque Indicating Ratchet Wrench (0-35 Ncm) comes w/ C9980	L-TIRW
Low Torque Indicating Ratchet Wrench Kit**	L-TIRWK
Standard ISO 1797 Adapter	C9980
Short ISO 1797 Adapter	C9981
ISO Latch, 4X4 Driver Adapter	ISO4X4A
Restorative Torque Tray	RTI2035TR



**L-TIRWK includes: L-TIRW, C9980, RASH3N, RASQ3N, RTI2035TR

Restorative Tooling Chart

External Hex

Square: torque 35 Ncm, Hexed: torque 20 Ncm. (May require square or slotted drivers)



Connection	Small Hex		Large Hex		Square		Slotted	
Length	Finger Driver	Driver Tip	Finger Driver	Driver Tip	Finger Driver	Driver Tip	Finger Driver	Driver Tip
17 mm	PHD00N		PHD02N	IHX1.2D*/ RASH3N	PSQD0N		PSD00	
24 mm	PHD01N	RASH2N	PHD03N	IHXL1.2D*/ RASH8N	PSQD1N	RASQ3N	PSD01	RASD1
30 mm		RASH7N				RASQ8N		RASD6
For use with:	Externally Hexed Cover Screws and Spline Cover Screws		Externally Hexed Encode Healing Abutments and Abutment Screws		Square Gold-Tite® Abutment Screws, Square Try-in Screws		Slotted Gold-Tite Screws, Core-Vent Comoponents, Slotted Implant Innovations Externally Hexed Restorative Components	

*GemLock®

OverdenSURE®

Torque 20 Ncm / Trilobe Connection

Certain® Implant System



Option

1



TWR (ODS-DRVR)

External Hex® Implant System



Option

2



L-TIRW (With ISO4X4A and ODS-DRVR)

For Full Arch Restorative Tooling, please see ZV2066.

For more information, visit ZimVie.com

ZimVie US Corp LLC
 4555 Riverside Drive
 Palm Beach Gardens, FL 33410
 Phone: +1-561-776-6700
 Fax: +1-561-776-1272
 dentalCS@ZimVie.com



Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. and its affiliates; and all products are manufactured by one or more of the dental subsidiaries of ZimVie US Corp LLC (Biomet 3i, LLC, etc.) and distributed by ZimVie US Corp LLC and its authorized marketing partners. Elos Torque Wrenches and Adapters are manufactured by Elos Medtech Pinol A/S. OverdenSURE and Azure tools are manufactured by Terrats Medical S.L. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie and is not for distribution in France. ZV3222-NAM REV A 01/26 ©2026 ZimVie. All rights reserved.

