


TSX® Immediate Molar Surgical Kit and Modules



- 1 **IMDR28L**
Pointed Starter Drill
2.8 mmD,
Long
- 2 **IMDR34TSX**
Tapered Drill
TSX 3.4 mmD
- 3 **IMDR38TSX**
Tapered Drill
TSX 3.8 mmD
- 4 **IMDR44TSX**
Tapered Drill
TSX 4.4 mmD
- 5 **IMDIR34**
Direction Indicator
3.4 mmD
- 6 **IMDIR38**
Direction Indicator
3.8 mmD
- 7 **IMDIR44**
Direction Indicator
4.4 mmD



- 8 **IMDR56**
Shaping Drill
5.4 mmD,
6 mmL
- 9 **IMDR66**
Shaping Drill
6 mmD,
6 mmL
- 10 **IMDR76**
Shaping Drill
7 mmD,
6 mmL
- 11 **IMDR86**
Shaping Drill
8 mmD,
6 mmL
- 12 **IMDR96**
Shaping Drill
9 mmD,
6 mmL
- 13 **IMDI56**
Depth Indicator
5.4 mmD,
6 mmL
- 14 **IMDI66**
Depth Indicator
6 mmD,
6 mmL
- 15 **IMDI76**
Depth Indicator
7 mmD,
6 mmL
- 16 **IMDI86**
Depth Indicator
8 mmD,
6 mmL
- 17 **IMDI96**
Depth Indicator
9 mmD,
6 mmL
- 18 **IMDR58**
Shaping Drill
5.4 mmD,
8 mmL
- 19 **IMDR68**
Shaping Drill
6 mmD,
8 mmL
- 20 **IMDR78**
Shaping Drill
7 mmD,
8 mmL
- 21 **IMDR88**
Shaping Drill
8 mmD,
8 mmL
- 22 **IMDR98**
Shaping Drill
9 mmD,
8 mmL
- 23 **IMDR510**
Shaping Drill
5.4 mmD,
10 mmL
- 24 **IMDR610**
Shaping Drill
6 mmD,
10 mmL
- 25 **IMDR710**
Shaping Drill
7 mmD,
10 mmL
- 26 **IMDR810**
Shaping Drill
8 mmD,
10 mmL
- 27 **IMDR910**
Shaping Drill
9 mmD,
10 mmL
- 28 **IMDR511**
Shaping Drill
5.4 mmD,
11.5 mmL
- 29 **IMDR611**
Shaping Drill
6 mmD,
11.5 mmL
- 30 **IMDR711**
Shaping Drill
7 mmD,
11.5 mmL
- 31 **IMDR811**
Shaping Drill
8 mmD,
11.5 mmL
- 32 **IMDR911**
Shaping Drill
9 mmD,
11.5 mmL
- 33 **IMDI5**
Depth Indicator
5.4 mmD,
8-11.5 mmL
- 34 **IMDI6**
Depth Indicator
6 mmD,
8-11.5 mmL
- 35 **IMDI7**
Depth Indicator
7 mmD,
8-11.5 mmL
- 36 **IMDI8**
Depth Indicator
8 mmD,
8-11.5 mmL
- 37 **IMDI9**
Depth Indicator
9 mmD,
8-11.5 mmL

- Tooling Module**
- 40 **RHX3.0DL**
Hex Driver, Latch-lock,
3.0 mm (D), 25 mm (L)
 - 41 **RHX3.0DS**
Hex Driver, Latch-lock,
3.0 mm (D), 21 mm (L)
 - 42 **RHX3.0RM**
Hex Driver, Ratchet,
3.0 mm (D), 22 mm (L)
 - 43 **RHX3.0RS**
Hex Driver, Ratchet,
3.0 mm (D), 17 mm (L)
 - 44 **HXLGR1.25**
GemLock
Hex Driver, Long
 - 45 **HXGR1.25**
GemLock
Hex Driver, Short
 - 46 **HX1.25D**
Latch-Lock
Hex Driver, Short,
1.25 mm(D), 23 mm(L)
 - 47 **Empty Grommets**
 - 48 **ZTIRW**
High Torque Indicating Ratchet (0-90 Ncm)

- Bone Tap Module (Optional)**
- 49 **ISOHEXADPT**
High Torque ISO
Hex Adapter
 - 50 **IMT76**
Bone Tap
7 mmD,
6 mmL
 - 51 **IMT86**
Bone Tap
8 mmD,
6 mmL
 - 52 **IMT96**
Bone Tap
9 mmD,
6 mmL
 - 53 **IMT7**
Bone Tap
7 mmD,
8-11.5 mmL
 - 54 **IMT8**
Bone Tap
8 mmD,
8-11.5 mmL
 - 55 **IMT9**
Bone Tap
9 mmD,
8-11.5 mmL

For more information, visit ZimVie.com



Images shown are representational in nature and may not be to scale or display the exact geometry or color of the actual components.

Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. ZTIRW is manufactured by Elos Medtech AB and distributed by ZimVie Dental. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZVINST0232 REV A 01/25 ©2025 ZimVie. All rights reserved.

