


Driva™ Drills Gold Series Usage Chart

TM, TSV® and Eztetic®

Use a dry erase marker to track the total number of osteotomy preparations that are performed with each drill. It is recommended that reusable drills be replaced after 15 uses. *Single use drills should be disposed of after each procedure.



*Single use drills must be disposed of after each procedure. This chart is intended to be used as a reference guide only. The number of uses per drill will vary and depends on a variety of factors including bone density encountered, proper handling and cleaning. Refer to Instructions for Use (IFU) for Biomet 3i Kits and Instruments (P-ZBDINSTRP) available at ifu.biomet3i.com for complete Cleaning and Sterilization Information.

NOTE: Surgical instruments are susceptible to damage for a variety of reasons including but not limited to: prolonged use, misuse, rough or improper handling, inadequate cleaning, etc. Care must be taken to avoid compromising the intended performance of the instrument.