

Dríva™ Gold Series Drills

Gold Metal Performance

Designed for high performance and excellent cutting efficiency, Dríva Gold Series Drills are reusable, externally irrigated, and compatible with existing surgical protocols for Trabecular Metal™, Tapered Screw-Vent® and Eztetic® Implants.



1 Efficient

- Deep cutting flutes designed for improved cutting efficiency compared to the original Dríva Drills
- Biocompatible ceramic coating¹ increases hardness and durability

2 Compatible

- Familiar soft- and dense-bone protocols are used for osteotomy preparation
- Black axial stripe on select drills indicates compatibility with Drill Stops

3 Convenient


- Solid construction facilitates easy cleaning and sterilization
- Gold color provides excellent visibility of grooves and bands that indicate drilling depth



ZIMMER BIOMET
Your progress. Our promise.®

Enhanced Cutting Efficiency¹

Original vs. Gold Series
(2.3 mm Driva Drill)



¹Data on file.

Ordering Information

Driva Gold Series Drills for Trabecular Metal and TSV® Implants

Diameter	Drill Length	
	16 mmL	22 mmL
2.3 mmD	TSV23SG	TSV23G
2.8 mmD	TSV28SG	TSV28G
3.4 mmD	TSV34SG	TSV34G
3.8 mmD	TSV38SG	TSV38G
5.1 mmD	TSV51SG	TSV51G

Driva Gold Series Step Drills for Trabecular Metal and TSV Implants

Diameter	Drill Length	
	16 mmL	22 mmL
3.4/2.8 mmD	TSV34D28SG	TSV34D28G
3.8/3.4 mmD	TSV38D34SG	TSV38D34G
4.4/3.8 mmD	TSV44D38SG	TSV44D38G
5.1/4.4 mmD	TSV51D44SG	TSV51D44G
5.7/5.1 mmD	TSV57D51SG	TSV57D51G

Driva Gold Series Surgical Kits for Trabecular Metal and TSV Implants

Description	Item No.
Gold Series Drill Set (New TSV Drills only, no tray)	TSVDRILLSG*
Tapered Screw-Vent Surgical Kit, Complete	TSVKITG
Tapered Screw-Vent Starter Kit (3.7 and 4.7 mmD only)	STRKITG
Surgical Tray and Driva Gold Series Drills for TSV only	TSVTRAYDG

*Available for a limited time only, while supplies last.

Driva Gold Series Module and Drills for Eztetic Implants

Description	Item No.
Eztetic Surgical Module, Driva Gold Series	NPMODG
Driva Gold Series Drill, 2.4 mmD, 22 mmL	EZT24G
Driva Gold Series Step Drill, 2.8/2.4 mmD, 22 mmL	EZT28D24G

Contact us at 1-800-342-5454 or visit zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB0874 REV A 09/19 ©2019 Zimmer Biomet. All rights reserved.

