

Ordering Information for Encode® Emergence Healing Abutments

Legacy 1.0 vs. Emergence



ORDERING INFORMATION

Encode Legacy 1.0 Sizes

Certain®

Restorative Platform	Emergence Profile	Collar Height			
		3 mmH	4 mmH	6 mmH	8 mmH
3.4 mmD	3.8 mmP	IEHA343	IEHA344	IEHA346	IEHA348
	5 mmP	IEHA353	IEHA354	IEHA356	IEHA358
4.1 mmD	4.1 mmP	IEHA443	IEHA444	IEHA446	IEHA448
	5 mmP	IEHA453	IEHA454	IEHA456	IEHA458
	6 mmP	IEHA463	IEHA464	IEHA466	IEHA468
5 mmD	5.6 mmP	IEHA553	IEHA554	IEHA556	IEHA558
	6 mmP	IEHA563	IEHA564	IEHA566	IEHA568
6 mmD	6.8 mmP	IEHA663	IEHA664	IEHA666	IEHA668



External Hex

Restorative Platform	Emergence Profile	Collar Height			
		3 mmH	4 mmH	6 mmH	8 mmH
3.4 mmD	3.8 mmP	EHA343	EHA344	EHA346	EHA348
	5 mmP	EHA353	EHA354	EHA356	EHA358
4.1 mmD	4.1 mmP	EHA443	EHA444	EHA446	EHA448
	5 mmP	EHA453	EHA454	EHA456	EHA458
	6 mmP	EHA463	EHA464	EHA466	EHA468
5 mmD	5.6 mmP	EHA553	EHA554	EHA556	EHA558
	6 mmP	EHA563	EHA564	EHA566	EHA568
6 mmD	6.8 mmP	EHA663	EHA664	EHA666	EHA668



TSV®

Restorative Platform	Emergence Profile	Collar Height		
		3 mmH	5 mmH	7 mmH
3.5 mmD	3.8 mmP	TEHA3383	TEHA3385	TEHA3387
	5 mmP	TEHA3503	TEHA3505	-
4.5 mmD	5 mmP	TEHA4503	TEHA4505	-
	5.6 mmP	TEHA4563	TEHA4565	TEHA4567
	6.0 mmP	TEHA4603	TEHA4605	-
5.7 mmD	6.8 mmP	TEHA5683	TEHA5685	-

Encode Emergence Sizes

Certain Connection Implants

Platform Diameter	Emergence Profile	Height			Quantity
		3	5	7	
3.4	3.8	IEEHA343	IEEHA345	IEEHA347	
	5.0	IEEHA353	IEEHA355	IEEHA357	
4.1	4.1	IEEHA443	IEEHA445	IEEHA447	
	5.0	IEEHA453	IEEHA455	IEEHA457	
	6.0	IEEHA463	IEEHA465	IEEHA467	
	7.5	IEEHA473	IEEHA475	IEEHA477	
5.0	5.0	IEEHA553	IEEHA555	IEEHA557	
	6.0	IEEHA563	IEEHA565	IEEHA567	
	7.5	IEEHA573	IEEHA575	IEEHA577	
6.0	6.0	IEEHA663	IEEHA665	-	
	6.8	IEEHA6683	IEEHA6685	-	
	7.5	IEEHA673	IEEHA675	-	

Tapered Screw-Vent® Connection Implants

Platform Diameter	Emergence Profile	Height			Quantity
		3	5	7	
3.5	3.8	TEEHA333	TEEHA335	TEEHA337	
	5.0	TEEHA353	TEEHA355	TEEHA357	
	6.5	TEEHA363	TEEHA365	TEEHA367	
4.5	4.5	TEEHA443	TEEHA445	TEEHA447	
	5.5	TEEHA453	TEEHA455	TEEHA457	
	6.5	TEEHA463	TEEHA465	TEEHA467	
	7.5	TEEHA473	TEEHA475	TEEHA477	
5.7	5.7	TEEHA553	TEEHA555	-	
	6.8	TEEHA563	TEEHA565	-	
	7.5	TEEHA573	TEEHA575	-	

References:

- Schubach P. The defense architecture of the human peri-implant mucosa: a histological study. J Prosthet Dent. 2007 Jun;97(6 Suppl):S15-25
- Abrahamsson. 1997 Journal of Clinical Periodontology 1997; 24: 565-572. The mucosal barrier following abutment dis/reconnection. An experimental study in dogs.

Ordering Information for Encode® Emergence Healing Abutments

Legacy 1.0 vs. Emergence



ORDERING INFORMATION

Healing Abutments

EP One-Piece

Restorative Platform	Emergence Profile	Collar Height				
		2 mmH	3 mmH	4 mmH	6 mmH	8 mmH
3.4 mmD	3.4 mmP	ISMHA32	ISMHA33	ISMHA34	ISMHA36	—
	3.8 mmP	IMHA32	IMHA33	IMHA34	IMHA36	—
	5 mmP	IMHA352	IMHA353	IMHA354	IMHA356	—
4.1 mmD	4.1 mmP	ISHA42	ISHA43	ISHA44	ISHA46	ISHA48
	5 mmP	ITHA52	ITHA53	ITHA54	ITHA56	ITHA58
	6 mmP	—	ITHA63	ITHA64	ITHA66	ITHA68
5.0 mmD	5 mmP	ISWHA52	ISWHA53	ISWHA54	ISWHA56	ISWHA58
	5.6 mmP	IWTH52	IWTH53	IWTH54	IWTH56	IWTH58
	6 mmP	IWTH562	IWTH563	IWTH564	IWTH566	IWTH568
6.0 mmD	6 mmP	ISWHA62	ISWHA63	ISWHA64	ISWHA66	ISWHA68
	6.8 mmP	IWTH62	IWTH63	IWTH64	IWTH66	IWTH68



Encode Emergence Sizes

Certain Connection Implants

Platform Diameter	Emergence Profile	Height			Quantity
		3	5	7	
3.4	3.8	IEEHA343	IEEHA345	IEEHA347	
	5.0	IEEHA353	IEEHA355	IEEHA357	
4.1	4.1	IEEHA443	IEEHA445	IEEHA447	
	5.0	IEEHA453	IEEHA455	IEEHA457	
	6.0	IEEHA463	IEEHA465	IEEHA467	
	7.5	IEEHA473	IEEHA475	IEEHA477	
5.0	5.0	IEEHA553	IEEHA555	IEEHA557	
	6.0	IEEHA563	IEEHA565	IEEHA567	
	7.5	IEEHA573	IEEHA575	IEEHA577	
6.0	6.0	IEEHA663	IEEHA665	-	
	6.8	IEEHA6683	IEEHA6685	-	
	7.5	IEEHA673	IEEHA675	-	

Healing Collars

Implant Platform	Emergence Profile	Cuff Height		
		3.0 mm	5.0 mm	7.0 mm
3.5 mmD	3.5 mmD (No Flare)	HC333	HC335	•
3.5 mmD	4.5 mmD	HC343	HC345	HC347
3.5 mmD	5.5 mmD	HC353	HC355	•
4.5 mmD	4.5 mmD (No Flare)	HC443	HC445	•
4.5 mmD	5.5 mmD	HC453	HC455	HC457
4.5 mmD	6.5 mmD	HC463	HC465	•
5.7 mmD	6.5 mmD	HC563	HC565	•



Tapered Screw-Vent® Connection Implants

Platform Diameter	Emergence Profile	Height			Quantity
		3	5	7	
3.5	3.8	TEEHA333	TEEHA335	TEEHA337	
	5.0	TEEHA353	TEEHA355	TEEHA357	
	6.5	TEEHA363	TEEHA365	TEEHA367	
4.5	4.5	TEEHA443	TEEHA445	TEEHA447	
	5.5	TEEHA453	TEEHA455	TEEHA457	
	6.5	TEEHA463	TEEHA465	TEEHA467	
	7.5	TEEHA473	TEEHA475	TEEHA477	
5.7	5.7	TEEHA553	TEEHA555	-	
	6.8	TEEHA563	TEEHA565	-	
	7.5	TEEHA573	TEEHA575	-	

For more information, visit ZimVie.com

ZimVie
 4555 Riverside Drive
 Palm Beach Gardens, FL 33410
 1-800-342-5454
 Phone: +1-561-776-6700
 Fax: +1-561-776-1272



Unless otherwise indicated, as referenced herein, all trademarks and intellectual property rights are the property of ZimVie Inc. or an affiliate; and all products are manufactured by one or more of the dental subsidiaries of ZimVie Inc. (Biomet 3i, LLC, Zimmer Dental, Inc., etc.) and marketed and distributed by ZimVie Dental and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. This material may not be copied or reprinted without the express written consent of ZimVie. ZV0702 REV A 02/23 ©2023 ZimVie. All rights reserved.

